


Quality Control with BeautifulJason



ISO 21561: Styrene-butadiene rubber (SBR) — Determination of the microstructure of solution- polymerized SBR

Area	Signal integration range
A	From 4,3 ppm to minimum intensity point around 5,0 ppm
B	From minimum intensity point around 5,0 ppm to minimum intensity point around 6,1 ppm
C	From minimum intensity point around 6,1 ppm to 7,7 ppm
TMS_{blank}	Integrated signal intensity of TMS in $CDCl_3$ containing TMS
CD_{blank}	From 6,1 ppm to 7,7 ppm in $CDCl_3$ containing TMS
TMS	Integrated signal intensity of TMS in S-SBR sample solution



$$C_{calib} = C - CD_{blank} \times (TMS/TMS_{blank}) \quad (1)$$

where

C_{calib} is the integrated signal intensity of area C compensated for the effect of $CHCl_3$ in $CDCl_3$.

3.6.3 Calculate the content of each microstructure component (trans and cis, and vinyl) of the butadiene portion and the styrene content, using Formulae (2) to (4):

$$S_m = \frac{(C_{calib}/5) \times 104}{(C_{calib}/5) \times 104 + (B/2 + A/4) \times 54} \times 100 \quad (2)$$

$$V = \frac{A/2}{B/2 + A/4} \times 100 \quad (3)$$

$$TC = \frac{B/2 - A/4}{B/2 + A/4} \times 100 \quad (4)$$

where

S_m is the styrene content of the S-SBR, in mass %;

V is the vinyl content of the butadiene portion of the S-SBR, in mol %;

TC is the trans and the cis content of the butadiene portion of the S-SBR, in mol %.

How to do this?

```
Windows PowerShell
Copyright (C) Microsoft Corporation. All rights reserved.


Install the latest PowerShell for new features and improvements! https://aka.ms/PSWindows

PS C:\Users\Botana.SERVDOMAIN2\Desktop\python_qc> python -m SBR_analysis_solvents subtraction
C Integral is :      21818.96
B Integral is :      12682.93
A Integral is :       3806.98
TMS Integral is :    13351.20
Solvent Integral is : 9805.93 at 7.26 ppm

=====
Styrene content is :          38.82 % mass
Vinyl content of butadiene is : 26.10 % molar
Trans and cis content of butadiene is : 73.90 % molar
=====

PS C:\Users\Botana.SERVDOMAIN2\Desktop\python_qc> |
```

Peak picking and creating rules



Script creation (I)

```
import beautifuljson as bjson

# python -m SBR_analysis_solventsubtraction

# Initialize variables
file1 = 'C:/Users/Botana.SERVDOMAIN2/Desktop/python_qc/SBR_proton-1-1.jdf'
rules_name = 'ISO_SBR'

# Initialize JASON object
jason = bjson.JASON()

# Open the file and apply the rules
with jason.create_document(file1, rules=rules_name) as doc:
    spec = doc.nmr_data[0]
```

Script creation (II)

```
# Extract the integrals and print the with 2 decimals
multiplet1 = spec.multiplets[0]
int_C_file1 = multiplet1.value_hz
print('C Integral is :      ', "{:.2f}".format(int_C_file1))
multiplet2 = spec.multiplets[1]
int_B_file1 = multiplet2.value_hz
print('B Integral is :      ', "{:.2f}".format(int_B_file1))
multiplet3 = spec.multiplets[2]
int_A_file1 = multiplet3.value_hz
print('A Integral is :      ', "{:.2f}".format(int_A_file1))
multiplet4 = spec.multiplets[3]
int_TMS_file1 = multiplet4.value_hz
print('TMS Integral is :    ', "{:.2f}".format(int_TMS_file1))
```

Script creation (III)

```
# Find the solvent peak and extract its area (assuming it is only one peak, otherwise keep adding the areas of each solvent peak and/or impurities)
solvent = None
for peak in spec.peaks:
    if peak.classification == bjson.NMRPeak.PeakClassification.NMRSolvent:
        solvent = peak
        break
int_solvent = solvent.area
print('Solvent Integral is : ', "{:.2f}".format(solvent.area), ' at ', "{:.2f}".format(solvent.pos[0]), ' ppm')

# Uncomment the below lines to save the document
# doc.close()
# doc.copy('C:/Users/Botana.SERVDOMAIN2/Desktop/python_qc/simple_integral_example.jjh5')
```

Script creation (IV)

```
# Subtract the solvent area from C integral
c_calib = int_C_file1 - int_solvent

# Calculate concentrations
s_m = c_calib*104/5*100/(c_calib*104/5+(int_B_file1/2+int_A_file1/4)*54)
v_ = int_A_file1/2*100/(int_B_file1/2+int_A_file1/4)
tc_ = (int_B_file1/2-int_A_file1/4)/(int_B_file1/2+int_A_file1/4)*100

print('')
print('=====')
print(' Styrene content is : ', "{:.2f}".format(s_m), ' % mass')
print(' Vinyl content of butadiene is : ', "{:.2f}".format(v_), ' % molar')
print(' Trans and cis content of butadiene is : ', "{:.2f}".format(tc_), ' % molar')
print('=====')
print('')
```