


STRUCTURE OF DELTA

Satoshi Sakurai
JEOL (UK) Ltd.

Basics

Operation PC and Console PC

We need to be aware of that both Operation PC and Console PC may store files to operate the system.






Global folder and Local folder

Basic rules of Delta software:

- In the case to use the same file name, the files in local folder (local files) have preference over the files which are stored in global folder (global files).
- The global files can not be directly edited.
- The local files can be customize independently for each Delta user.

Main components

Main components involved in the operation of the system are



Sample conditions

Experimental file


Process list

Solvent information

Automation script

Delta V4

Delta V4 (V3)



Delta V4



Operation PC

Global folder
Local folder

Experimental file

Automation script

Process list

Solvent information




Console PC

Sample conditions

Solvent information

Delta V4 (V3)



Delta V4



Operation PC




Experimental file


Automation script

Process list

Solvent information



Delta V4 (V3)



Delta V4

```
管理者: コマンド プロンプト
eclipse
file_digest.sha1sum
files
gamma_def
global_op_experiment_tool_ex.jpj
global_op_parameter_database.jpj
global_op_peak_analyzer.jpj
global_op_plot_param_ex.jpj
global_op_pulse_width_analyzer.jpj
global_op_sequence_file_generator.jpj
global_op_super_automation_queue.jpj
global_op_super_automation_tool.jpj
global_operator_library.jpj
machine.alarms
machine.announce
machine.authorization
machine.config
machine.config_0
machine.config_02
machine.config_orig
machine.connect
machine.defaults
machine.defaults65531
machine.defaults65532
machine.defaults65533
machine.defaults65534
machine.defaults65535
machine.dm_shims
machine.idle_lock
phantom:/eclipse#
machine.post_exp65533
machine.post_exp65534
machine.post_exp65535
machine.pre_exp
machine.pre_exp65531
machine.pre_exp65532
machine.pre_exp65533
machine.pre_exp65534
machine.pre_exp65535
machine.rr_shims
machine.wiring
parameter.config
probes
queue.hours
queue.priorities
queue.sample_id
shapes
solvent.def
solvent.def_0
solvent.def_02
solvent_reference.def
start
start.org
stop
stop.org
templates
usr
verify_install
```




Console PC

Sample conditions

Solvent information

By Command prompt, it is available to modify...

Delta V4 (V3)



Delta V4



Operation PC

Global folder
Local folder

Experimental file

Automation script

Process list

Solvent information




Console PC

Sample conditions

Solvent information

Delta V5

Delta V5



Delta V5



Operation PC



Console PC

Global folder
Local folder

Experimental file

Automation script

Process list

Global folder
Local folder

Experimental file


Sample conditions

Automation script

Solvent information

Process list

The merit to use local folder of operation PC



Delta V5



Operation PC

Global folder
Local folder

Experimental file

Automation script

Process list

- It is specific for computer.
- It is possible to restrict to Window's account.

The note to use local folder of V5 for operation PC

There are three different “local” folder!!

- “local” 1 (C:¥Users¥”account”¥Documents¥JEOL¥Delta 5.2)

i.e. in the case of “Delta” account,

→ C:¥Users¥Delta¥Documents¥JEOL¥Delta 5.2

Experimental file → C:¥Users¥”account name”¥Documents¥JEOL¥Delta 5.2¥experiments

Automation script → C:¥Users¥”account name”¥Documents¥JEOL¥Delta 5.2¥automation

Process list → C:¥Users¥”account name”¥Documents¥JEOL¥Delta 5.2¥process_lists

It is for personal use.

- “local” 2 (C:¥Users¥”account”¥AppData¥Local¥JEOL¥Delta 5.2)

*AppDta folder is hidden folder.

It is a folder to store the temporal setting file.

For example, when you change the preferences, system.jnv file is stored here.

- “local” 3 (C:¥Program Files¥Common Files¥JEOL¥Delta 5.2)

This folder is restricted for administrator account to modify.

It is available to restrict the some setting by making local “.cfg” file here.

→explain more by next slide.


- Ref. Global folder (C:¥Program Files¥JEOL¥Delta 5.2.0.app¥Contents¥Resources¥global)




Operation PC

How does it affect?

i.e. to restrict the language to "Russian".



The merit to use local folder of console PC




Delta V5



Operation PC



2nd operation PC



Console PC

Global folder
Local folder

- Experimental file
- Automation script
- Process list
- Sample conditions
- Solvent information

- You can use the same files from other computer.
- You can easily share the original files with colleagues.

The note to use local folder of V5 for console PC

There are two different “local” folder!!

- “local” 1 (C:¥Program Files¥Common Files¥JEOL¥Control 5.2)

Experimental file → C:¥Program Files¥Common Files¥JEOL¥Control 5.2¥experiments

Automation script → C:¥Program Files¥Common Files¥JEOL¥Control 5.2¥automation

Process list → C:¥Program Files¥Common Files¥JEOL¥Control 5.2¥process_lists

Sample conditions, Solvent information → C:¥Program Files¥Common Files¥JEOL¥Control 5.2¥configuration

It is possible to access by any account



- “local” 2 (C:¥Program Files¥Common Files¥JEOL¥Control 5.2¥users¥”account”)

Experimental file, Automation script, Process list, Solvent information

→ C:¥Program Files¥Common Files¥JEOL¥Control 5.2¥users¥”account”

*It is not valid for sample conditions.

i.e. in the case of “Delta” account,

→ C:¥Program Files¥Common Files¥JEOL¥Control 5.2¥users¥Delta

It is specific for Delta account



- Ref. Global folder (C:¥Program Files¥JEOL¥Control 5.2.0.app¥Contents¥Resources¥global)



Console PC

How to make local files for Delta V5 (For console PC)



How to make a local file

Experimental file

Automation script


Process list

Upload via file browser of Delta software



Edit via preference tool of Delta software


Sample conditions



Console PC

Edit via Remote Maintenance Tool

Solvent information



Upload via file browser of Delta software


- “local” 1 (C:¥Program Files¥Common Files¥JEOL¥Control 5.2)

Experimental file

Automation script

Process list

i.e. upload an experimental file.



Upload via file browser of Delta software


- “local” 2 (C:¥Program Files¥Common Files¥JEOL¥Control 5.2¥users¥”account”)

Experimental file


Automation script

Process list


i.e. upload an experimental file for Satoshi account.




1. Click “Open file and choose tool”.




3. Select [File]–[Upload File...].




5. Do authentication by “Satoshi”.



2. Select target file.



4. Click “Connect As...” button.




6. Select “Satoshi” folder.
7. Click “OK” button.

How to access to the local file

i.e. access to experimental file.


From “local” 1

On “Open Experiment”, select “Global” folder from left box, and develop “all files” folder. There is uploaded experimental file with description of “Local” there.



From “local” 2

On “Open Experiment”, select “Local” folder from left box. There is uploaded experimental file there.



How does it affect?



Delta



Satoshi

“local” 1 folder

“local” 2 folder (Satoshi)

none

none

Proton Carbon COSY DEPT

Available Methods

- Standard
 - COSY
 - Carbon
 - DEPT
 - HMBC
 - HMOC
 - HSQC
 - NOESY
 - Proton
 - ROESY
 - TOCSY
 - edited DEPT

Proton Carbon COSY DEPT

Available Methods

- Standard
 - COSY
 - Carbon
 - DEPT
 - HMBC
 - HMOC
 - HSQC
 - NOESY
 - Proton
 - ROESY
 - TOCSY
 - edited DEPT

```

AUTOMATION VERSION 2;+
PURPOSE "A collection of common Methods for basic experimentation.;"
INCLUDE "Utilities" TO DOMAIN util;+
-----
METHOD Local1 ( OUT highres_proton : DATA ) IS+
  PURPOSE "Routine Proton single pulse experiment";+
  EXPERIMENT Proton IS+
    SAVE AS "${SAMPLE}_Proton";+
    COLLECT "proton";+
  END EXPERIMENT;+
END METHOD Local1;+
LEOF
    
```

none

Proton Carbon COSY DEPT

Available Methods

- Standard
 - Local1

Proton Carbon COSY DEPT

Available Methods

- Standard
 - Local1

```

AUTOMATION VERSION 2;+
PURPOSE "A collection of common Methods for basic experimentation.;"
INCLUDE "Utilities" TO DOMAIN util;+
-----
METHOD Local1 ( OUT highres_proton : DATA ) IS+
  PURPOSE "Routine Proton single pulse experiment";+
  EXPERIMENT Proton IS+
    SAVE AS "${SAMPLE}_Proton";+
    COLLECT "proton";+
  END EXPERIMENT;+
END METHOD Local1;+
LEOF
    
```

```

AUTOMATION VERSION 2;+
PURPOSE "A collection of common Methods for basic experimentation.;"
INCLUDE "Utilities" TO DOMAIN util;+
-----
METHOD Local2 ( OUT highres_proton : DATA ) IS+
  PURPOSE "Routine Proton single pulse experiment";+
  EXPERIMENT Proton IS+
    SAVE AS "${SAMPLE}_Proton";+
    COLLECT "proton";+
  END EXPERIMENT;+
END METHOD Local2;+
LEOF
    
```

Proton Carbon COSY DEPT

Available Methods

- Standard
 - Local1

Proton Carbon COSY DEPT

Available Methods



- Standard
 - Local2

Edit via preference tool of Delta software


- “local” 1 (C:¥Program Files¥Common Files¥JEOL¥Control 5.2)

Sample conditions

1. Login to spectrometer by “Console” account.
2. Open the “Hardware Configuration” from [Config]-[Instrument Preferences].
3. On “Hardware Configuration” window, open “Edit Preference” from [Parameters]-[Edit...].



4. Put the commands as next slide.



Edit via preference tool of Delta software


- “local” 1 (C:¥Program Files¥Common Files¥JEOL¥Control 5.2)

Sample conditions


Originally, sample conditions are defined on “[sample_fields.jnd](#)” which is in C:¥Program Files¥JEOL¥Control 5.2.0.app¥Contents¥Resources¥global¥configuration.

```
GROUP.GROUP(WIDGET) = "INPUT_BOX"  
GROUP.GROUP(CONTAINER_TYPE_INCLUDE) = "NUMBER"  
  
-----  
- LIQUIDS  
-----  
  
LIQUIDS.TEMP_SET(PRETTY)+PARAMETER.TEMP_SET = "temp_set" —"degrees"  
LIQUIDS.TEMP_SET(DEFAULT) = 25[dC]  
LIQUIDS.TEMP_SET(WIDGET) = "CONSTRAINED_INPUT_BOX"  
LIQUIDS.TEMP_SET(ASYNCHRONOUS) = FALSE  
LIQUIDS.TEMP_SET(FILTER) = ""  
LIQUIDS.TEMP_SET(AUTOMATIC) = TRUE  
  
LIQUIDS."SPIN_SET"(PRETTY)+PARAMETER.SPIN_SET = "spin_set" —"spin rate"  
LIQUIDS."SPIN_SET"(DEFAULT) = 15[Hz]  
LIQUIDS."SPIN_SET"(WIDGET) = "CONSTRAINED_INPUT_BOX"  
LIQUIDS."SPIN_SET"(ASYNCHRONOUS) = FALSE  
LIQUIDS."SPIN_SET"(FILTER) = ""  
LIQUIDS."SPIN_SET"(AUTOMATIC) = TRUE  
  
LIQUIDS.SPIN_STATE(PRETTY)+PARAMETER.SPIN_STATE = "spin_state" —"spinner"  
LIQUIDS.SPIN_STATE(DEFAULT) = "SPIN_ON"  
LIQUIDS.SPIN_STATE(WIDGET) = "SELECT_ENUMERATION"  
LIQUIDS.SPIN_STATE(ASYNCHRONOUS) = FALSE  
LIQUIDS.SPIN_STATE(FILTER) = ""  
LIQUIDS.SPIN_STATE(AUTOMATIC) = TRUE  
-LIQUIDS.SPIN_STATE(PRETTY) = "spinner"  
-LIQUIDS.SPIN_STATE(TYPE) = "BOOLEAN"  
-LIQUIDS.SPIN_STATE(DEFAULT) = TRUE  
-LIQUIDS.SPIN_STATE(WIDGET) = "TOGGLE_BUTTON"  
  
LIQUIDS.TEMP_STATE(PRETTY)+PARAMETER.TEMP_STATE = "temp_state"  
- temperature  
LIQUIDS.TEMP_STATE(DEFAULT) = "TEMP_OFF"  
LIQUIDS.TEMP_STATE(WIDGET) = "SELECT_ENUMERATION"  
LIQUIDS.TEMP_STATE(ASYNCHRONOUS) = FALSE  
LIQUIDS.TEMP_STATE(FILTER) = ""  
LIQUIDS.TEMP_STATE(AUTOMATIC) = TRUE
```

If I want to set the temperature control with 30 [dC] for all liquid sample as default.



It requires rebooting of Delta software of control.



How does it affect?

The new "sample_fields.jnv" file was saved in "local"1 folder.

```
— sample_fields  
— C:\Program Files\Common Files\JEOL\Control 5.2\configuration\sample_fields.jnv  
liquids.spin_state(default) = "SPIN OFF"  
liquids.temp_set(default) = 30[dC]  
liquids.temp_state(default) = "TEMP ON"  
standard."load shims"(automatic) = TRUE  
standard."load shims"(default) = "System"  
—#Signature 7BF7D06C4A40AFD5402E80CDC79260AA80BE785FD JEOL 10-OCT-2017  
17:41:23 GMT 螟乘函髯
```

Before

No.	Sample Name	Solvent	Slot	Kind	Shared	Verified	Err
1	test	DMSO-D6	1	Liquids	<input type="checkbox"/>	<input checked="" type="checkbox"/>	
2		NONE	0	Liquids	<input type="checkbox"/>	<input type="checkbox"/>	

Attributes	Value
spin_set	15[Hz]
spin_state	SPIN OFF
temp_set	25[dC]
temp_state	TEMP OFF

After



No.	Sample Name	Solvent	Slot	Kind	Shared	Verified	Err
1	test	DMSO-D6	1	Liquids	<input type="checkbox"/>	<input checked="" type="checkbox"/>	
2		NONE	0	Liquids	<input type="checkbox"/>	<input type="checkbox"/>	

Attributes	Value
spin_set	15[Hz]
spin_state	SPIN OFF
temp_set	30[dC]
temp_state	TEMP ON

Edit via Remote Maintenance Tool

Solvent information


1. Open the Remote Maintenance Tool (Instrument).
2. Login by console privilege user.
3. Move to “Solvent Management” session.
4. Put in the all information about solvent, like a “Solvent name” and etc.
5. Click the “Add” button.



Remote Maintenance Tool

Structure of pulse program


Script of Experimental file



```
-----  
--  
-- Experiment Source Code --  
-- Delta NMR Experiment & Machine Control Interface --  
--  
-- Copyright (c) 2000 JEOL Ltd --  
-- All Rights Reserved --  
-----  
--  
-- HELP.eng: absorption pn_hsqc using club with X-decoupling  
-- Category: basic, 2D, hsqc, liquids  
-- File name : hsqc  
--  
-- Sequence name :absorption pn_hsqc using club with X-decoupling  
--  
-- Reference :  
-- H.Hu,A.J.Shaka, J.Magn.Reson.,136,54-629(1999)  
--  
-- V.A.Mandelsham,H.Hu,A.J.Shaka,Magn.Reson.Chem.,36,S17-S28(1998)  
--  
-- R.R. Ernst, G. Bodenhausen, A. Wokaun "Principles of Nuclear Magnetic  
-- Resonance in One and Two Dimensions" Clarendon Press,  
-- Oxford University Press, 1987.  
--  
-- G.E. Martin, A.S. Zektzer, "Two Dimensional NMR Methods for Establishing  
-- NMR Connectivity, VCH, 1988.  
--  
-- J.K.M Sanders, B Hunter "Modern NMR Spectroscopy, A Guide for Chemists"  
-- Oxford University Press, 1987.  
--  
-- Parameters  
-- x_pulse : 90[deg] obs pulse width  
-- x_atn : attenuator of x_pulse  
-- y_pulse : 90[deg] irr pulse width  
-- y_atn : attenuator of y_pulse  
--  
-- x_fm_pulse : 180 deg FM pulse width along x axis  
-- x_fm_atn : attenuator of x_fm_pulse  
--  
-- y_fm_pulse : 180 deg FM pulse width along y axis  
-- y_fm_atn : attenuator of y_fm_pulse
```

script

Script of Experimental file



```
comment 8 =? " *** Pulse Field Gradient ***";
grad_selection =? if y_domain = "Carbon13" then "13C = 2.5:1.5:1" else
if y_domain = "Nitrogen15" then "15N = 5.94:3.94:1" else
if y_domain = "Silicon29" then "29Si = 3.52:1.52:1" else
if y_domain = "Phosphorus31" then "31P = 2.24:0.24:1" else
"none";

include "services";
x_gamma = _get_gamma( x_domain );
y_gamma = _get_gamma( y_domain );
temp1 = x_gamma / y_gamma;

grad_1 => 1[ms], help "duration of grad_1";
grad_1_amp => 0.075[T/m], -0.3[T/m] -> 0.3[T/m]:0.0003[T/m], help "amplitude of
grad_1";
grad_1a_amp = grad_1_amp/2;
grad_1b_amp = -1*grad_1_amp/2;

grad_2 => grad_1, help "duration of grad_2";
grad_2_amp => 0.045[T/m], -0.3[T/m] -> 0.3[T/m]:0.0003[T/m], help "amplitude of
grad_2";
grad_2a_amp = -1*grad_2_amp/2;
grad_2b_amp = grad_2_amp/2;

grad_3 => grad_1, help "duration of grad_3";
grad_3_amp => (grad_1_amp+grad_2_amp)/temp1, help "amplitude of grad_3";

grad_shape => "SINE", fg_shape_names, help "gradient shape";
grad_recover => 0.1[ms], help "gradient recovery time";

-y_p1_correction =? 360 * ((x_fm_pulse + 1.28*y_pulse+40[us])* y_sweep);

include "y_dec_param_carbon13";
include "obs_dante_presat";
include "obs_presat_time";
include "tri_presat_homo_dec";
include "presat_relaxation_delay_calc";
include "pulse";

phase_x = {0};
phase_y = {90};
phase_qd = {2(0, 180)}.ystep(180%2);
```

script

.jxi files

By command of “include”, it downloads “.jxi” file from global file (C:\Program Files\JEOL\Control 5.2.0.app\Contents\Resources\global\experiments).

```
—y_p1_correction =? 360 * ((x_fm_pulse + 1.28*y_pulse+40[us])* y_sweep);  
  
include "y_dec_param_carbon13";  
include "obs_dante_presat";  
include "obs_presat_time";  
include "tri_presat_homo_dec";  
include "presat_relaxation_delay_calc";  
include "pulse";
```




```
—y_dec_param_carbon13.jxi  
—  
comment_102 =? "*** im_decoupling ***";  
  
im_decoupling => TRUE, help "select TRUE or FALSE for decoupling or non_decoupling";  
  
when im_decoupling do  
  im_noise => "MPF8" ("CW" "DIPS12" "WAUGH" "WALTZ" "GARP"  
    "MPF5" "MPF6" "MPF7" "MPF8" "MPF9" "MPF10"  
    "WURST_30" "WURST_40" "WURST_50" "WURST_70"), help "decoupler  
  noise modulation",  
end when;  
  
include "y_dec_param";  
|
```

y_dec_param_carbon13.jxi

It is necessary to put the .jxi file into the same folder with experimental file in order to include it into experimental file.

Develop the include file into experimental file

When you want to develop the parameters which are defined in include file into experimental file, it is possible by following process.



```
grad_3_amp      => (grad_1_amp + grad_2_amp) / temp1, help "amplitude of grad_
3";
grad_shape      => "SINE", fg_shape_names, help "gradient shape";
grad_recover    => 0.1[ms], help "gradient recovery time";
comment_102     =? "*** ir_decoupling ***";
ir_decoupling   => TRUE, help "select TRUE or FALSE for decoupling or
non_decoupling";
when ir_decoupling do
  ir_noise      => "MPF8", ("CW" "DIPS12" "WAUGH" "WALTZ" "GARF"
"MPF5" "MPF6" "MPF7" "MPF8" "MPF9" "MPF10" "WURST_30" "WURST_40"
"WURST_50" "WURST_70"), help "decoupler noise modulation";
end when;
ir_domain       = y_domain;
when ir_decoupling do
  ir_wurst      = if ir_noise = "WURST_30"
then TRUE
else if ir_noise = "WURST_40"
then TRUE
else if ir_noise = "WURST_50"
then TRUE
else if ir_noise = "WURST_70"
then TRUE
else FALSE, help "Select TRUE or FALSE for adiabatic
decoupling";
end when;
when ir_wurst do
  ir_wurst_b1_default =? if y_domain = "Fluorine19"
then 5[kHz]
else 2.5[kHz];
ir_wurst_b1      => ir_wurst_b1_default;
ir_wurst_pwidth  = 1 / (4 * ir_wurst_b1);
ir_decoupling_factor_value = if x_domain = "Carbon13"
then if ir_domain = "Fluorine19"
then 0.7
else 0.9
else 0.9;
```


Summary


Structure of operation PC




Operation PC



Console PC



Structure of console PC



Summary

By understanding the file folder structure of the and its priority, it is possible to customize the environment for each account.

Thank you for your kind attention

