



DOSY automatic setup

Adolfo Botana, PhD
JEOL UK
12th September 2022

DOSY Discovery Day 2022

Self-Diffusion and molecular weight (SEGWE)

- Approximation for the correlation between molecular weight and self-diffusion coefficient
- More advanced than previous estimations such as $r_H \propto \sqrt[3]{MW}$
- Takes into account the molecular interactions of solvent and solutes
- **Calculator available, input: temperature, solvent and MW or D**

k_B : Boltzmann constant

T : temperature

η : viscosity

ρ_{eff} : effective density of the molecule
(packing effects, geometry,
solvation and flexibility)

MW : molecular weight of the
molecule

MW_S : molecular weight of the solvent

N_A : Avogadro number

$$D = \frac{k_B T \left(\frac{3 \sqrt[3]{\frac{MW_S}{MW}}}{2} + \frac{1}{1 + \sqrt[3]{\frac{MW_S}{MW}}} \right)}{6\pi\eta \sqrt[3]{\frac{3MW}{4\pi\rho_{\text{eff}} N_A}}}$$

R. Evans et al, *Angew. Chem. Int. Ed.* (2013), 52: 3199–3202

R. Evans et al, *Anal. Chem.* (2018), 90 (6), 3987–3994

R. Evans, *Prog Nucl Magn Reson Spectrosc* (2020), 117, 33-69

<https://www.nmr.chemistry.manchester.ac.uk/?q=node/432>

Calculating gradient length (Stejskal-Tanner equation)

$$S = S_0 e^{-D\gamma^2\delta^2G^2\Delta'}$$

$$\delta = \sqrt{-\frac{\ln S/S_0}{D\gamma^2G^2\Delta'}}$$

S : signal amplitude [0.1 S_0]

S_0 : signal amplitude without diffusion [S_0]

D : diffusion coefficient [SEGWE]

γ : gyromagnetic ratio (dependent on x_domain) [γ_H]

δ : gradient pulse width (**delta** or **smallDelta**)

G : gradient amplitude (**g**) [maximum usable gradient strength]

Δ' : corrected diffusion time

Δ : diffusion time (**diffusion_time**) [0.1 s]

DOSY steps (automatic_setup)

1. Automatic setting of chosen sample parameters (temperature, temperature delay, shimming, tuning*,...)
2. Automatic estimation of diffusion coefficient using SEGWE with the **solute MW**
3. Automatic determination of gradient length to achieve 90% signal attenuation (0.1 s diffusion time and maximum usable gradient strength)

Info

Starting Job 'DOSY'
Estimated Diffusion coefficient is 10.98386*10E-10 m²/s
Delta set to 1.8041[ms] to achieve 90% signal attenuation

4. Automatic determination of 90 degrees pulse and setup (**optional**)
5. Automatic determination of sample convection (**optional**), by evaluating if a delay imbalance of 30ms leads to a signal loss of more than 20%.

Info Post-experiment Default Initialization
There is convection in this sample (signal ratio of 0.59402. Proceeding with convection compensated experiment
Digital Filter Factor is 12

Info Post-experiment Default Initialization
Convection test indicates there is no significant convection in this sample (signal ratio of 0.98005)
Digital Filter Factor is 12

6. Automatic acquisition of DOSY with/without convection compensation
7. Automatic DOSY plot automatically generation

dosy_scans	16
dosy_relaxation_delay	3[s]
solute_MW	300
calculate_proton_90	<input type="checkbox"/>
convection_check	<input checked="" type="checkbox"/>
automatic_setup	<input checked="" type="checkbox"/>
diffusion_time	0.1[s]
smallDelta	2[ms]
predefined_array	<input checked="" type="checkbox"/>
g	5[T/m], 0.28993[T/m], 0.3[T/m]}
array_type	Linear ↕
g_max	0.3[T/m]
g_min	0.03[T/m]
g_points	16
log_base	2

DOSY steps (semi-automatic setup)

1. Automatic setting of chosen sample parameters (temperature, temperature delay, shimming, tuning*,...)
2. Choose **diffusion time** and **gradient length**
3. If predefined_array, type in/copy-paste **gradient strengths**
If not predefined, setup array as linear / exponential / logarithmic with the weakest PFG strength, strongest PFG strength, (logarithm base if logarithmic array) and number of PFG values
3. Automatic determination of 90 degrees pulse and setup (**optional**)
4. Automatic determination of sample convection (**optional**), by evaluating if a delay imbalance of 30ms leads to a signal loss of more than 20%.


Info Post-experiment Default Initialization
There is convection in this sample (signal ratio of 0.59402. Proceeding with convection compensated experiment
Digital Filter Factor is 12

Info Post-experiment Default Initialization
Convection test indicates there is no significant convection in this sample (signal ratio of 0.98005)
Digital Filter Factor is 12

5. Automatic acquisition of DOSY with/out convection compensation
6. Automatic DOSY plot automatically generation

▶	dosy_scans	16
▶	dosy_relaxation_delay	3[s]
▶	solute_MW	300
▶	calculate_proton_90	<input type="checkbox"/>
▶	convection_check	<input checked="" type="checkbox"/>
▶	automatic_setup	<input type="checkbox"/>
▶	diffusion_time	0.1[s]
▶	smallDelta	2[ms]
▶	predefined_array	<input type="checkbox"/>
▶	g	5[T/m], 0.28993[T/m], 0.3[T/m]}
▶	array_type	Linear ↓
▶	g_max	0.3[T/m]
▶	g_min	0.03[T/m]
▶	g_points	16
▶	log_base	2

Results from automatic_setup (Ethyl-indanone)



PDF automatically generated after clicking the Submit Job button with the DOSY automation method

Note PDF printout is generated using multivariate analysis, results can be misleading

Further analysis with other methods is recommended, particularly for more complex samples

Results from automatic_setup (Quinine + Geraniol + Camphene)

MW set to 324

Info Reached Lock State 'AUTOLOCK'
Estimated Diffusion coefficient is $7.78834 \times 10^{-10} \text{ m}^2/\text{s}$
Delta set to 2.14247[ms] to achieve 90% signal attenuation



1H 90 degree pulse width for QCG estimated as 6.12549[us]

Room temperature

Convection test indicates there is no significant convection in this sample (signal ratio of 0.9784)

25C

There is convection in this sample (signal ratio of 0.35024).
Proceeding with convection compensated experiment



Results from automatic_setup (Quinine + Geraniol + Camphene)

MW set to 136

Info

Reached Lock State 'AUTOLOCK'
 Estimated Diffusion coefficient is $11.71103 \times 10^{-10} \text{ m}^2/\text{s}$
 Delta set to 1.74719[ms] to achieve 90% signal attenuation

