

---

# Permanent PDA Configuration with LCMS-9050

## 1. Overview

This guide outlines the setup and maintenance of an alternative configuration where the **column outlet is permanently connected through the SPD-M30A PDA**, and the **PDA outlet is routed to the switching valve**. The switching valve is used to select or deselect the LCMS-9050.

---

## 2. System Configuration

### 2.1 Plumbing Setup

- Connect the **column outlet** directly to the **inlet of the SPD-M30A PDA**.
- Connect the **PDA outlet** to the **switching valve**, which then routes flow either to the **LCMS-9050** or to waste, depending on the selected mode.

### 2.2 Instrument Configuration

- When running in **MS-Only mode**, ensure:
    - The **valve is set** to bypass the PDA output to the LCMS-9050.
    - The **PDA is turned off** and not included in the '9050-Only' instrument configuration.
    - This renders the PDA a **passive flow-through element**.
- 

## 3. Advantages of This Configuration

- **Simplified Tandem Operation:** Enables easy transition to a tandem PDA-MS configuration (referred to as the "Forbidden" setup).
  - **Improved Flow-Cell Maintenance:** Continuous solvent flow through the PDA reduces the risk of:
    1. **Contamination**
    2. **Degradation from stop-flow with D2 lamp on**
    3. **Drying out during storage**
- 

## 4. Potential Disadvantages

- **Post-Column Dead Volume:** Adds ~10–15 µL, which may slightly reduce chromatographic resolution when using UPLC columns. However, in DIA mode, this may be beneficial due to slight peak broadening.
- **Particulate Contamination:** Risk of blockage in the flow-cell. Use a **particulate filter** upstream of the PDA.
- **Background Signal:** Running PDA at high concentrations may increase background in MS-only runs due to residual deposits.
- **Method Hygiene:** Strictly segregate methods by configuration to avoid errors or damage.

---

## 5. Maintenance Guidelines

### 5.1 Flow-Cell Care

- **Avoid Drying Out:** Always store the cell filled with **100% methanol**.
- **Monthly Flushing:** Flush the cell with pure methanol once per month to prevent degradation.
- **Avoid Stop-Flow with Lamp On:** Ensure the **D2 lamp is off** when there is no flow. If the pump is on and the PDA is accidentally powered, the shutter may open and degrade the cell.


### 5.2 Operational Best Practices

- **Keep PDA Off in MS-Only Runs:** Prevents unnecessary wear and background signal.
  - **Use MS-Compatible Solvents:** Ensures compatibility and longevity of both PDA and MS systems.
  - **Monitor for Noise:** Interference patterns (450–500 nm) may appear due to refractive index changes—this is normal and related to the capillary's internal reflection.
- 

## 6. Notes on Flow-Cell Design

- The SPD-M30A uses an **85 mm capillary** with **Total Internal Reflection** to transmit light.
  - The flow-cell is a **consumable**; internal components are not user-replaceable.
  - Most failures are due to **drying out** or **burnt contamination** from improper use.
- 


### Alternative LCMS-9050 or HPLC-PDA Flow Selection Configuration



Note that the Left-hand Configuration and Valve Position would now remain the same when using the SPD-M30A and LCMS-9050, in Tandem.

The difference is which *Instrument Configuration Profile* is opened from LabSolutions, and whether the SPD-M30A is powered On (for PDA + 9050 use), or Off (for LCMS-9050-Only)

## LCMS-9050 or HPLC-PDA Flow Selection



Connect outlet from HPLC column to *Port 1*, on the valve (inlet *from* HPLC column).

Select *Switching-Valve Position '1'* in the software (currently the Default for LCMS-9050 Configuration)

Apply flow to column, and determine from which port ('2', or '6') solvent flows

Stop flow, and connect this port (shown here as '2' – but depends on rotor orientation) to the tubing leading to the LCMS-9050 *Divert Valve*

Download Valve Position '0' to the column oven, from Realtime, and confirm solvent flows from other port (shown here as '6')

Connect this port to PDA

Load 'LCMS-9050' Instrument Configuration. Check and, if necessary, **change All Methods used with this Config.** to Valve Position '1'

Load 'HPLC-PDA' Instrument Configuration. Check and, if necessary, **change All Methods used with this Config.** to Valve Position '0'

**Clearly label methods, as to which instrument configuration they are intended for, and only use methods for the correct configuration**